

TCP11L2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00177-1

Capture Antibody Information

Catalog Number: 83227-1-PBS	Clone ID: 240096F6	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: t-complex 11 (mouse)-like 2
Isotype: IgG	GenBank: BC033617	Gene ID: 255394
Purification Method: Protein A purification	Immunogen Catalog Number: Ag11185	

Detection Antibody Information

Catalog Number: 83227-2-PBS	Clone ID: 240096H2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: t-complex 11 (mouse)-like 2
Isotype: IgG	GenBank: BC033617	Gene ID: 255394
Purification Method: Protein A purification	Immunogen Catalog Number: Ag11185	

Applications

Tested Applications: Cytometric bead array	Range: 1.25-20 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP00177-1 targets TCP11L2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: TCP11L2 Recombinant antibody, PBS Only (Capture) 83227-1-PBS (240096F6). 100 µg. Concentration 1 mg/mL.

Detection antibody: TCP11L2 Recombinant antibody, PBS Only (Detector) 83227-2-PBS (240096H2). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve ofMP00177-1, TCP11L2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83227-1-PBS. Detection antibody: 83227-2-PBS. Standard: Ag11185. Range: 1.25-20 ng/mL