

CEP350-Specific Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00178-2**

Capture Antibody Information

Catalog Number:
83208-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230320D4

Reactivity:
Human

GenBank:
NM_014810

Conjugate:
Unconjugated

Full name:
centrosomal protein 350kDa

Gene ID:
9857

Detection Antibody Information

Catalog Number:
83208-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230320A5

Reactivity:
Human

GenBank:
NM_014810

Conjugate:
Unconjugated

Full name:
centrosomal protein 350kDa

Gene ID:
9857

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00178-2 targets CEP350-Specific in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CEP350-Specific Recombinant antibody, PBS Only (Capture/Detector) 83208-2-PBS (230320D4). 100 µg. Concentration 1 mg/mL.

Detection antibody: CEP350-Specific Recombinant antibody, PBS Only (Detector) 83208-1-PBS (230320A5). 100 µg. Concentration 1 mg/mL.

Alternative CEP350-Specific matched antibody pairs: MP00178-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve ofMP00178-2, CEP350-Specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83208-2-PBS. Detection antibody: 83208-1-PBS. Standard: SY00642. Range: 1.25-20 ng/mL