

# E-cadherin Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00187-1

## Capture Antibody Information

<b>Catalog Number:</b> 80541-5-PBS	<b>Clone ID:</b> 240159C10	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> Human	<b>Full name:</b> cadherin 1, type 1, E-cadherin (epithelial)
<b>Isotype:</b> IgG	<b>GenBank:</b> BC141838	<b>Gene ID:</b> 999
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag14973	

## Detection Antibody Information

<b>Catalog Number:</b> 80541-2-PBS	<b>Clone ID:</b> 240159A5	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> Human	<b>Full name:</b> cadherin 1, type 1, E-cadherin (epithelial)
<b>Isotype:</b> IgG	<b>GenBank:</b> BC141838	<b>Gene ID:</b> 999
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag14973	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.625-20 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
------------------------------------------------------	---------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------

## Product Information

MP00187-1 targets E-cadherin in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: E-cadherin Recombinant antibody, PBS Only (Capture) 80541-5-PBS (240159C10). 100 µg. Concentration 1 mg/mL.

Detection antibody: E-cadherin Recombinant antibody, PBS Only (Detector) 80541-2-PBS (240159A5). 100 µg. Concentration 1 mg/mL.

Alternative E-cadherin matched antibody pairs: MP00187-2, MP00187-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP00187-1, E-cadherin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 80541-5-PBS. Detection antibody: 80541-2-PBS. Standard: Ag14973. Range: 0.625-20 ng/mL.