

NSD3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00189-2

Capture Antibody Information

Catalog Number: 83116-2-PBS	Clone ID: 230499C2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: Wolf-Hirschhorn syndrome candidate 1-like 1
Isotype: IgG	GenBank: BC012059	Gene ID: 54904
Purification Method: Protein A purification	Immunogen Catalog Number: Ag1894	

Detection Antibody Information

Catalog Number: 83116-4-PBS	Clone ID: 230499F12	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: Wolf-Hirschhorn syndrome candidate 1-like 1
Isotype: IgG	GenBank: BC012059	Gene ID: 54904
Purification Method: Protein A purification	Immunogen Catalog Number: Ag1894	

Applications

Tested Applications: Cytometric bead array	Range: 0.625-20 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00189-2 targets NSD3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NSD3 Recombinant antibody, PBS Only (Capture) 83116-2-PBS (230499C2). 100 µg. Concentration 1 mg/mL.

Detection antibody: WHSC1L1 Recombinant antibody, PBS Only (Detector) 83116-4-PBS (230499F12). 100 µg. Concentration 1 mg/mL.

Alternative NSD3 matched antibody pairs: MP00189-1, MP00189-3, MP00189-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00189-2, WHSC1L1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83116-2-PBS. Detection antibody: 83116-4-PBS. Standard: Ag1894. Range: 0.625-20 ng/mL