

SLC6A9 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00243-1**

Capture Antibody Information

Catalog Number: 83150-3-PBS	Clone ID: 240010H2	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: solute carrier family 6 (neurotransmitter transporter, glycine), member 9
Isotype: IgG	GenBank: NM_201649	Gene ID: 6536
Purification Method: Protein A purification	Immunogen Catalog Number: Ag34437	

Detection Antibody Information

Catalog Number: 83150-1-PBS	Clone ID: 240010A7	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: solute carrier family 6 (neurotransmitter transporter, glycine), member 9
Isotype: IgG	GenBank: NM_201649	Gene ID: 6536
Purification Method: Protein A purification	Immunogen Catalog Number: Ag34437	

Applications

Tested Applications: Cytometric bead array	Range: 0.625-20 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00243-1 targets SLC6A9 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.
Capture antibody: SLC6A9 Recombinant antibody, PBS Only (Capture) 83150-3-PBS (240010H2). 100 µg. Concentration 1 mg/mL.

Detection antibody: SLC6A9 Recombinant antibody, PBS Only (Detector) 83150-1-PBS (240010A7). 100 µg. Concentration 1 mg/mL.

Alternative SLC6A9 matched antibody pairs: MP00243-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00243-1, SLC6A9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83150-3-PBS. Detection antibody: 83150-1-PBS. Standard: Ag34437. Range: 0.625-20 ng/mL