

SLC6A9 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00243-2**

Capture Antibody Information

Catalog Number:
83150-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240010E8
Reactivity:
Human
GenBank:
NM_201649
Immunogen Catalog Number:
Ag34437

Conjugate:
Unconjugated
Full name:
solute carrier family 6
(neurotransmitter transporter,
glycine), member 9
Gene ID:
6536

Detection Antibody Information

Catalog Number:
83150-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240010A7
Reactivity:
Human
GenBank:
NM_201649
Immunogen Catalog Number:
Ag34437

Conjugate:
Unconjugated
Full name:
solute carrier family 6
(neurotransmitter transporter,
glycine), member 9
Gene ID:
6536

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00243-2 targets SLC6A9 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC6A9 Recombinant antibody, PBS Only (Capture) 83150-2-PBS (240010E8). 100 μ g. Concentration 1 mg/mL.

Detection antibody: SLC6A9 Recombinant antibody, PBS Only (Detector) 83150-1-PBS (240010A7). 100 μ g. Concentration 1 mg/mL.

Alternative SLC6A9 matched antibody pairs: MP00243-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00243-2, SLC6A9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83150-2-PBS. Detection antibody: 83150-1-PBS. Standard: Ag34437. Range: 1.25-20 ng/mL.