

TIMP1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00245-1**

Capture Antibody Information

Catalog Number:
83151-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230329E8
Reactivity:
Human
GenBank:
BC000866

Conjugate:
Unconjugated
Full name:
TIMP1 metalloproteinase inhibitor 1
Gene ID:
7076
Affinity:
 $K_D = 2.16 \times 10^{-11}$
 $K_{off} = 4.92 \times 10^{-5}$
 $K_{on} = 2.28 \times 10^6$

Detection Antibody Information

Catalog Number:
83151-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230329F11
Reactivity:
Human
GenBank:
BC000866

Conjugate:
Unconjugated
Full name:
TIMP1 metalloproteinase inhibitor 1
Gene ID:
7076
Affinity:
 $K_D = 1.75 \times 10^{-11}$
 $K_{off} = 5.52 \times 10^{-6}$
 $K_{on} = 3.15 \times 10^5$

Applications

Tested Applications:
Cytometric bead array

Range:
0.625–20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00245-1 targets TIMP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TIMP1 Recombinant antibody, PBS Only (Capture) 83151-3-PBS (230329E8). 100 µg. Concentration 1 mg/mL.

Detection antibody: TIMP1 Recombinant antibody, PBS Only (Detector) 83151-2-PBS (230329F11). 100 µg. Concentration 1 mg/mL.

Alternative TIMP1 matched antibody pairs: MP00245-2, MP00245-3, MP50084-1, MP50084-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00245-1, TIMP1 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 83151-3-PBS. Detection antibody: 83151-2-PBS. Standard: Eg0638. Range: 0.625-20 ng/mL