

RANK/TNFRSF11A Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00249-1

Capture Antibody Information

Catalog Number:
83253-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240106E7

Reactivity:
human

GenBank:
NM_003839

Immunogen Catalog Number:
Ag29425

Conjugate:
Unconjugated

Full name:
tumor necrosis factor receptor superfamily, member 11a, NFKB activator

Gene ID:
8792

Detection Antibody Information

Catalog Number:
83253-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240106B10

Reactivity:
human

GenBank:
NM_003839

Immunogen Catalog Number:
Ag29425

Conjugate:
Unconjugated

Full name:
tumor necrosis factor receptor superfamily, member 11a, NFKB activator

Gene ID:
8792

Applications

Tested Applications:
Cytometric bead array

Range:
2.5-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00249-1 targets RANK/TNFRSF11A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RANK Recombinant antibody, PBS Only (Capture) 83253-2-PBS (240106E7). 100 µg. Concentration 1 mg/mL.

Detection antibody: RANK Recombinant antibody, PBS Only (Detector) 83253-1-PBS (240106B10). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00249-1, TNFRSF 11A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83253-2-PBS. Detection antibody: 83253-1-PBS. Standard: Ag29425. Range: 2.5-40 ng/mL.