

PPAR Gamma Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00256-1

Capture Antibody Information

Catalog Number:
81490-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230374C2
Reactivity:
human
GenBank:
BC006811
Immunogen Catalog Number:
Ag10005

Conjugate:
Unconjugated
Full name:
peroxisome proliferator-activated
receptor gamma
Gene ID:
5468

Detection Antibody Information

Catalog Number:
81490-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230374E7
Reactivity:
human
GenBank:
BC006811
Immunogen Catalog Number:
Ag10005

Conjugate:
Unconjugated
Full name:
peroxisome proliferator-activated
receptor gamma
Gene ID:
5468

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-20 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP00256-1 targets PPAR Gamma in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PPAR Gamma Recombinant antibody, PBS Only (Capture) 81490-3-PBS (230374C2). 100 µg.
Concentration 1 mg/mL.

Detection antibody: PPAR Gamma Recombinant antibody, PBS Only (Detector) 81490-4-PBS (230374E7). 100 µg.
Concentration 1 mg/mL.

Alternative PPAR Gamma matched antibody pairs: MP00256-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00256-1, PPAR Gamma Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81490-3-PBS. Detection antibody: 81490-4-PBS. Standard: Ag10005. Range: 1.25-20 ng/mL