

EFHD2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00269-1

Capture Antibody Information

Catalog Number:
83264-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240185C3
Reactivity:
Human
GenBank:
BC023611
Immunogen Catalog Number:
Ag35084

Conjugate:
Unconjugated
Full name:
EF-hand domain family, member D2
Gene ID:
79180

Detection Antibody Information

Catalog Number:
83264-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240185G3
Reactivity:
Human
GenBank:
BC023611
Immunogen Catalog Number:
Ag35084

Conjugate:
Unconjugated
Full name:
EF-hand domain family, member D2
Gene ID:
79180

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00269-1 targets EFHD2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: EFHD2 Recombinant antibody, PBS Only (Capture/Detector) 83264-1-PBS (240185C3). 100 µg. Concentration 1 mg/mL.

Detection antibody: EFHD2 Recombinant antibody, PBS Only (Detector) 83264-4-PBS (240185G3). 100 µg. Concentration 1 mg/mL.

Alternative EFHD2 matched antibody pairs: MP00269-2, MP00269-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00269-1, EFHD2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83264-1-PBS. Detection antibody: 83264-4-PBS. Standard: Ag35084. Range: 1.25-40 ng/mL.