

# SLC6A12 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00289-3**

## Capture Antibody Information

**Catalog Number:**  
83195-4-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
240021D9

**Reactivity:**  
human

**GenBank:**  
BC126215

**Immunogen Catalog Number:**  
Ag18778

**Conjugate:**  
Unconjugated

**Full name:**  
solute carrier family 6 (neurotransmitter transporter, betaine/GABA), member 12

**Gene ID:**  
6539

## Detection Antibody Information

**Catalog Number:**  
83195-1-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
240021A8

**Reactivity:**  
human

**GenBank:**  
BC126215

**Immunogen Catalog Number:**  
Ag18778

**Conjugate:**  
Unconjugated

**Full name:**  
solute carrier family 6 (neurotransmitter transporter, betaine/GABA), member 12

**Gene ID:**  
6539

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-20 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00289-3 targets SLC6A12 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC6A12 Recombinant antibody, PBS Only (Capture) 83195-4-PBS (240021D9). 100 µg. Concentration 1 mg/mL.

Detection antibody: SLC6A12 Recombinant antibody, PBS Only (Detector) 83195-1-PBS (240021A8). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


## Storage

**Storage:**  
Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00289-3, SLC6A12 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83195-4-PBS. Detection antibody: 83195-1-PBS. Standard: Ag18778. Range: 0.625-20 ng/mL.