

PROCA1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00294-1

Capture Antibody Information

Catalog Number:
83193-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230538C10
Reactivity:
Human
GenBank:
BC029574
Immunogen Catalog Number:
Ag10851

Conjugate:
Unconjugated
Full name:
proline-rich cyclin A1-interacting protein
Gene ID:
147011

Detection Antibody Information

Catalog Number:
83193-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
230538F2
Reactivity:
Human
GenBank:
BC029574
Immunogen Catalog Number:
Ag10851

Conjugate:
Unconjugated
Full name:
proline-rich cyclin A1-interacting protein
Gene ID:
147011

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00294-1 targets PROCA1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PROCA1 Recombinant antibody, PBS Only (Capture) 83193-1-PBS (230538C10). 100 µg. Concentration 1 mg/mL.

Detection antibody: PROCA1 Recombinant antibody, PBS Only (Detector) 83193-2-PBS (230538F2). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve ofMP00294-1, PROCA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83193-1-PBS. Detection antibody: 83193-2-PBS. Standard: Ag10851. Range: 1.25-40 ng/mL