

CDX2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00304-1**

Capture Antibody Information

Catalog Number:

83125-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

230289G1

Reactivity:

Human

GenBank:

NM_001265

Conjugate:

Unconjugated

Full name:

caudal type homeobox 2

Gene ID:

1045

Detection Antibody Information

Catalog Number:

83125-4-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

230289G12

Reactivity:

Human

GenBank:

NM_001265

Conjugate:

Unconjugated

Full name:

caudal type homeobox 2

Gene ID:

1045

Applications

Tested Applications:

Cytometric bead array

Range:

0.313-10 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00304-1 targets CDX2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CDX2 Recombinant antibody, PBS Only (Capture) 83125-1-PBS (230289G1). 100 µg. Concentration 1 mg/mL.

Detection antibody: CDX2 Recombinant antibody, PBS Only (Detector) 83125-4-PBS (230289G12). 100 µg. Concentration 1 mg/mL.

Alternative CDX2 matched antibody pairs: MP00304-2, MP50121-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:

Store at -80°C.

Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve ofMP00304-1, CDX2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83125-1-PBS. Detection antibody: 83125-4-PBS. Standard: SY01336. Range: 0.313-10 ng/mL