

GAPDHS Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00323-1

Capture Antibody Information

| | | |
|---|--|--|
| Catalog Number: 83290-1-PBS | Clone ID: 240277E4 | Conjugate: Unconjugated |
| Host: Rabbit | Reactivity: human | Full name: glyceraldehyde-3-phosphate dehydrogenase, spermatogenic |
| Isotype: IgG | GenBank: BC036373 | Gene ID: 26330 |
| Purification Method: Protein A purification | Immunogen Catalog Number: Ag4988 | |

Detection Antibody Information

| | | |
|---|--|--|
| Catalog Number: 83290-2-PBS | Clone ID: 240277H3 | Conjugate: Unconjugated |
| Host: Rabbit | Reactivity: human | Full name: glyceraldehyde-3-phosphate dehydrogenase, spermatogenic |
| Isotype: IgG | GenBank: BC036373 | Gene ID: 26330 |
| Purification Method: Protein A purification | Immunogen Catalog Number: Ag4988 | |

Applications

| | | |
|--|---|---|
| Tested Applications: Cytometric bead array | Range: 3.13-100 ng/mL (Cytometric Bead Array) | Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results. |
|--|---|---|

Product Information

MP00323-1 targets GAPDHS in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GAPDHS Recombinant antibody, PBS Only (Capture) 83290-1-PBS (240277E4). 100 µg. Concentration 1 mg/mL.

Detection antibody: GAPDHS Recombinant antibody, PBS Only (Detector) 83290-2-PBS (240277H3). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00323-1, GAPDHS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83290-1-PBS. Detection antibody: 83290-2-PBS. Standard: Ag4988. Range: 3.13-100 ng/mL