

# DDX47 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00325-1

## Capture Antibody Information

**Catalog Number:**  
83287-2-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
240193E2

**Reactivity:**  
human

**GenBank:**  
BC068009

**Immunogen Catalog Number:**  
Ag11953

**Conjugate:**  
Unconjugated

**Full name:**  
DEAD (Asp-Glu-Ala-Asp) box polypeptide 47

**Gene ID:**  
51202

## Detection Antibody Information

**Catalog Number:**  
83287-1-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
240193A6

**Reactivity:**  
human

**GenBank:**  
BC068009

**Immunogen Catalog Number:**  
Ag11953

**Conjugate:**  
Unconjugated

**Full name:**  
DEAD (Asp-Glu-Ala-Asp) box polypeptide 47

**Gene ID:**  
51202

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-80 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00325-1 targets DDX47 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DDX47 Recombinant antibody, PBS Only (Capture) 83287-2-PBS (240193E2). 100 µg. Concentration 1 mg/mL.

Detection antibody: DDX47 Recombinant antibody, PBS Only (Detector) 83287-1-PBS (240193A6). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00325-1, DDX47 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83287-2-PBS. Detection antibody: 83287-1-PBS. Standard: Ag11953. Range: 0.625-80 ng/mL