

SMARCA4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00334-3**

Capture Antibody Information

Catalog Number:
83310-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240172H1

Reactivity:
Human

GenBank:
BC150298

Immunogen Catalog Number:
Ag16256

Conjugate:
Unconjugated

Full name:
SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 4

Gene ID:
6597

Detection Antibody Information

Catalog Number:
83310-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240172B10

Reactivity:
Human

GenBank:
BC150298

Immunogen Catalog Number:
Ag16256

Conjugate:
Unconjugated

Full name:
SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 4

Gene ID:
6597

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00334-3 targets SMARCA4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SMARCA4 Recombinant antibody, PBS Only (Capture) 83310-4-PBS (240172H1). 100 µg. Concentration 1 mg/mL.

Detection antibody: SMARCA4 Recombinant antibody, PBS Only (Detector) 83310-3-PBS (240172B10). 100 µg. Concentration 1 mg/mL.

Alternative SMARCA4 matched antibody pairs: MP00334-1, MP00334-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00334-3, SMARCA4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83310-4-PBS. Detection antibody: 83310-3-PBS. Standard: Ag16256. Range: 0.625-40 ng/mL.