

TRIM34 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00340-3**

Capture Antibody Information

Catalog Number:
83316-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240180E2
Reactivity:
Human
GenBank:
NM_001003827
Immunogen Catalog Number:
Ag34444

Conjugate:
Unconjugated
Full name:
tripartite motif-containing 34
Gene ID:
53840

Detection Antibody Information

Catalog Number:
83316-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240180D10
Reactivity:
Human
GenBank:
NM_001003827
Immunogen Catalog Number:
Ag34444

Conjugate:
Unconjugated
Full name:
tripartite motif-containing 34
Gene ID:
53840

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00340-3 targets TRIM34 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TRIM34 Recombinant antibody, PBS Only (Capture) 83316-2-PBS (240180E2). 100 µg. Concentration 1 mg/mL.

Detection antibody: TRIM34 Recombinant antibody, PBS Only (Detector) 83316-3-PBS (240180D10). 100 µg. Concentration 1 mg/mL.

Alternative TRIM34 matched antibody pairs: MP00340-1, MP00340-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00340-3, TRIM34 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83316-2-PBS. Detection antibody: 83316-3-PBS. Standard: Ag34444. Range: 1.25-40 ng/mL.