

ADRB2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00341-1

Capture Antibody Information

Catalog Number:
83317-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240186G7

Reactivity:
Human

GenBank:
BC012481

Immunogen Catalog Number:
Ag3715

Conjugate:
Unconjugated

Full name:
adrenergic, beta-2-, receptor, surface

Gene ID:
154

Detection Antibody Information

Catalog Number:
83317-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240186D12

Reactivity:
Human

GenBank:
BC012481

Immunogen Catalog Number:
Ag3715

Conjugate:
Unconjugated

Full name:
adrenergic, beta-2-, receptor, surface

Gene ID:
154

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00341-1 targets ADRB2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ADRB2 Recombinant antibody, PBS Only (Capture) 83317-2-PBS (240186G7). 100 µg. Concentration 1 mg/mL.

Detection antibody: ADRB2 Recombinant antibody, PBS Only (Detector) 83317-3-PBS (240186D12). 100 µg. Concentration 1 mg/mL.

Alternative ADRB2 matched antibody pairs: MP00341-2, MP00341-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00341-1, ADRB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83317-2-PBS. Detection antibody: 83317-3-PBS. Standard: Ag3715. Range: 0.156-20 ng/mL