

# ZO-1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00349-3**

## Capture Antibody Information

<b>Catalog Number:</b> 82870-2-PBS	<b>Clone ID:</b> 240150A1	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> human	<b>Full name:</b> tight junction protein 1 (zona occludens 1)
<b>Isotype:</b> IgG	<b>GenBank:</b> BC111712	<b>Gene ID:</b> 7082
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag33182	

## Detection Antibody Information

<b>Catalog Number:</b> 82870-4-PBS	<b>Clone ID:</b> 240150B7	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> human	<b>Full name:</b> tight junction protein 1 (zona occludens 1)
<b>Isotype:</b> IgG	<b>GenBank:</b> BC111712	<b>Gene ID:</b> 7082
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag33182	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 2.5-80 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

## Product Information

MP00349-3 targets ZO-1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ZO1 Recombinant antibody, PBS Only (Capture) 82870-2-PBS (240150A1). 100 µg. Concentration 1 mg/mL.

Detection antibody: ZO1 Recombinant antibody, PBS Only (Detector) 82870-4-PBS (240150B7). 100 µg. Concentration 1 mg/mL.

Alternative ZO-1 matched antibody pairs: MP00349-1, MP00349-2, MP00349-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00349-3, ZO-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82870-2-PBS. Detection antibody: 82870-4-PBS. Standard: Ag33182. Range: 2.5-80 ng/mL.