

# SLC39A10 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00367-1

## Capture Antibody Information

**Catalog Number:**  
83358-3-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
240261G2

**Reactivity:**  
Human

**GenBank:**  
BC101516

**Immunogen Catalog Number:**  
Ag20099

**Conjugate:**  
Unconjugated

**Full name:**  
solute carrier family 39 (zinc transporter), member 10

**Gene ID:**  
57181

## Detection Antibody Information

**Catalog Number:**  
83358-2-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
240261G3

**Reactivity:**  
Human

**GenBank:**  
BC101516

**Immunogen Catalog Number:**  
Ag20099

**Conjugate:**  
Unconjugated

**Full name:**  
solute carrier family 39 (zinc transporter), member 10

**Gene ID:**  
57181

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-80ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00367-1 targets SLC39A10 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC39A10 Recombinant antibody, PBS Only (Capture) 83358-3-PBS (240261G2). 100 µg. Concentration 1 mg/mL.

Detection antibody: SLC39A10 Recombinant antibody, PBS Only (Detector) 83358-2-PBS (240261G3). 100 µg. Concentration 1 mg/mL.

Alternative SLC39A10 matched antibody pairs: MP00367-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.

**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00367-1, SLC39A10 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83358-3-PBS. Detection antibody: 83358-2-PBS. Standard: Ag20099. Range: 0.625-80ng/mL.