

# KIAA1199 / CEMIP Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00373-1**

## Capture Antibody Information

<b>Catalog Number:</b> 83373-3-PBS	<b>Clone ID:</b> 240372C3	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> human	<b>Full name:</b> KIAA1199 / CEMIP
<b>Isotype:</b> IgG	<b>GenBank:</b> BC020256	<b>Gene ID:</b> 57214
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag15527	

## Detection Antibody Information

<b>Catalog Number:</b> 83373-1-PBS	<b>Clone ID:</b> 240372B7	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> human	<b>Full name:</b> KIAA1199 / CEMIP
<b>Isotype:</b> IgG	<b>GenBank:</b> BC020256	<b>Gene ID:</b> 57214
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag15527	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 0.625-80 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

## Product Information

MP00373-1 targets KIAA1199 / CEMIP in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: KIAA1199 Recombinant antibody, PBS Only (Capture) 83373-3-PBS (240372C3). 100 µg. Concentration 1 mg/mL.

Detection antibody: KIAA1199 Recombinant antibody, PBS Only (Capture/Detector) 83373-1-PBS (240372B7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00373-1, KIAA1199 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83373-3-PBS. Detection antibody: 83373-1-PBS. Standard: Ag15527. Range: 0.625-80 ng/mL.