

# RBM20 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00379-2**

## Capture Antibody Information

**Catalog Number:**  
83363-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240345C5  
**Reactivity:**  
human  
**GenBank:**  
NM\_001134363  
**Immunogen Catalog Number:**  
Ag34322

**Conjugate:**  
Unconjugated  
**Full name:**  
RNA binding motif protein 20  
**Gene ID:**  
282996

## Detection Antibody Information

**Catalog Number:**  
83363-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240345B7  
**Reactivity:**  
human  
**GenBank:**  
NM\_001134363  
**Immunogen Catalog Number:**  
Ag34322

**Conjugate:**  
Unconjugated  
**Full name:**  
RNA binding motif protein 20  
**Gene ID:**  
282996

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
1.25-80 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00379-2 targets RBM20 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RBM20 Recombinant antibody, PBS Only (Capture) 83363-3-PBS (240345C5). 100 µg. Concentration 1 mg/mL.

Detection antibody: RBM20 Recombinant antibody, PBS Only (Capture/Detector) 83363-1-PBS (240345B7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00379-2, RBM20 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83363-3-PBS. Detection antibody: 83363-1-PBS. Standard: Ag34322. Range: 1.25-80 ng/mL.