

RBM20 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00379-4**

Capture Antibody Information

Catalog Number: 83363-2-PBS	Clone ID: 240345F1	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: RNA binding motif protein 20
Isotype: IgG	GenBank: NM_001134363	Gene ID: 282996
Purification Method: Protein A purification	Immunogen Catalog Number: Ag34322	

Detection Antibody Information

Catalog Number: 83363-1-PBS	Clone ID: 240345B7	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: RNA binding motif protein 20
Isotype: IgG	GenBank: NM_001134363	Gene ID: 282996
Purification Method: Protein A purification	Immunogen Catalog Number: Ag34322	

Applications

Tested Applications: Sandwich ELISA	Range: 78.1-5000 pg/mL (Sandwich ELISA)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	---

Product Information

MP00379-4 targets RBM20 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: RBM20 Recombinant antibody, PBS Only (Capture/Detector) 83363-2-PBS (240345F1). 100 µg. Concentration 1 mg/mL.

Detection antibody: RBM20 Recombinant antibody, PBS Only (Capture/Detector) 83363-1-PBS (240345B7). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Sandwich ELISA standard curve of MP00379-4, Human RBM20 Recombinant Matched Antibody Pair - PBS only. 83363-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag34322. 83363-1-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.



The mean RBM20 concentration was determined to be 1.8 ng/mL in A549 cell extract based on a 2.0 mg/mL extract load.