

CRIP1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00382-2

Capture Antibody Information

Catalog Number:
83340-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240176C5

Reactivity:
Human

GenBank:
BC002738

Immunogen Catalog Number:
Ag7588

Conjugate:
Unconjugated

Full name:
cysteine-rich protein 1 (intestinal)

Gene ID:
1396

Detection Antibody Information

Catalog Number:
83340-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240176G2

Reactivity:
Human

GenBank:
BC002738

Immunogen Catalog Number:
Ag7588

Conjugate:
Unconjugated

Full name:
cysteine-rich protein 1 (intestinal)

Gene ID:
1396

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00382-2 targets CRIP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CRIP1 Recombinant antibody, PBS Only (Capture) 83340-1-PBS (240176C5). 100 μ g. Concentration 1 mg/mL.

Detection antibody: CRIP1 Recombinant antibody, PBS Only (Capture/Detector) 83340-3-PBS (240176G2). 100 μ g. Concentration 1 mg/mL.

Alternative CRIP1 matched antibody pairs: MP00382-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00382-2, CRIP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83340-1-PBS. Detection antibody: 83340-3-PBS. Standard: Ag7588. Range: 0.156-20 ng/mL