

C1QTNF2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00385-1

Capture Antibody Information

Catalog Number:
83343-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240330E5

Reactivity:
Human

GenBank:
BC011699

Immunogen Catalog Number:
Ag9183

Conjugate:
Unconjugated

Full name:
C1q and tumor necrosis factor related protein 2

Gene ID:
114898

Detection Antibody Information

Catalog Number:
83343-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240330F8

Reactivity:
Human

GenBank:
BC011699

Immunogen Catalog Number:
Ag9183

Conjugate:
Unconjugated

Full name:
C1q and tumor necrosis factor related protein 2

Gene ID:
114898

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00385-1 targets C1QTNF2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: C1QTNF2 Recombinant antibody, PBS Only (Capture) 83343-3-PBS (240330E5). 100 µg. Concentration 1 mg/mL.

Detection antibody: C1QTNF2 Recombinant antibody, PBS Only (Detector) 83343-2-PBS (240330F8). 100 µg. Concentration 1 mg/mL.

Alternative C1QTNF2 matched antibody pairs: MP00385-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00385-1, C1QTNF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83343-3-PBS. Detection antibody: 83343-2-PBS. Standard: Ag9183. Range: 0.313-40 ng/mL