

MTF1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00386-1

Capture Antibody Information

Catalog Number:
83344-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240336C12

Reactivity:
Human

GenBank:
BC014454

Immunogen Catalog Number:
Ag21848

Conjugate:
Unconjugated

Full name:
metal-regulatory transcription factor 1

Gene ID:
4520

Detection Antibody Information

Catalog Number:
83344-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240336F7

Reactivity:
Human

GenBank:
BC014454

Immunogen Catalog Number:
Ag21848

Conjugate:
Unconjugated

Full name:
metal-regulatory transcription factor 1

Gene ID:
4520

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00386-1 targets MTF1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MTF1 Recombinant antibody, PBS Only (Capture/Detector) 83344-1-PBS (240336C12). 100 µg. Concentration 1 mg/mL.

Detection antibody: MTF1 Recombinant antibody, PBS Only (Detector) 83344-3-PBS (240336F7). 100 µg. Concentration 1 mg/mL.

Alternative MTF1 matched antibody pairs: MP00386-2, MP00386-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00386-1, MTF1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83344-1-PBS. Detection antibody: 83344-3-PBS. Standard: Ag21848. Range: 0.156-20 ng/mL