

LYPD3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00387-1

Capture Antibody Information

Catalog Number: 83345-3-PBS	Clone ID: 240347D8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: LY6/PLAUR domain containing 3
Isotype: IgG	GenBank: NM_014400	Gene ID: 27076
Purification Method: Protein A purification	Immunogen Catalog Number: Ag34961	

Detection Antibody Information

Catalog Number: 83345-2-PBS	Clone ID: 240347F4	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: LY6/PLAUR domain containing 3
Isotype: IgG	GenBank: NM_014400	Gene ID: 27076
Purification Method: Protein A purification	Immunogen Catalog Number: Ag34961	

Applications

Tested Applications: Cytometric bead array	Range: 0.156-20 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00387-1 targets LYPD3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: LYPD3 Recombinant antibody, PBS Only (Capture/Detector) 83345-3-PBS (240347D8). 100 µg. Concentration 1 mg/mL.

Detection antibody: LYPD3 Recombinant antibody, PBS Only (Detector) 83345-2-PBS (240347F4). 100 µg. Concentration 1 mg/mL.

Alternative LYPD3 matched antibody pairs: MP00387-2, MP00387-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00387-1, LYPD3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83345-3-PBS. Detection antibody: 83345-2-PBS. Standard: Ag34961. Range: 0.156-20 ng/mL