

DUSP6 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00391-1

Capture Antibody Information

Catalog Number:
83355-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240303A7
Reactivity:
Human
GenBank:
BC003562
Immunogen Catalog Number:
Ag0705

Conjugate:
Unconjugated
Full name:
dual specificity phosphatase 6
Gene ID:
1848

Detection Antibody Information

Catalog Number:
83355-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240303C4
Reactivity:
Human
GenBank:
BC003562
Immunogen Catalog Number:
Ag0705

Conjugate:
Unconjugated
Full name:
dual specificity phosphatase 6
Gene ID:
1848

Applications

Tested Applications:
Cytometric bead array

Range:
2.5-80ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00391-1 targets DUSP6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DUSP6 Recombinant antibody, PBS Only (Capture) 83355-1-PBS (240303A7). 100 µg. Concentration 1 mg/mL.

Detection antibody: DUSP6 Recombinant antibody, PBS Only (Detector) 83355-2-PBS (240303C4). 100 µg. Concentration 1 mg/mL.

Alternative DUSP6 matched antibody pairs: MP00391-2, MP00391-3, MP00391-4, MP00391-5, MP00391-6

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00391-1, DUSP6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83355-1-PBS. Detection antibody: 83355-2-PBS. Standard: Ag0705. Range: 2.5-80ng/mL.