

# PCARE Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00392-1

## Capture Antibody Information

**Catalog Number:**  
83364-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240352B6  
**Reactivity:**  
Human  
**GenBank:**  
NM\_001029883  
**Immunogen Catalog Number:**  
Ag34774

**Conjugate:**  
Unconjugated  
**Full name:**  
chromosome 2 open reading frame 71  
**Gene ID:**  
388939

## Detection Antibody Information

**Catalog Number:**  
83364-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240352A7  
**Reactivity:**  
Human  
**GenBank:**  
NM\_001029883  
**Immunogen Catalog Number:**  
Ag34774

**Conjugate:**  
Unconjugated  
**Full name:**  
chromosome 2 open reading frame 71  
**Gene ID:**  
388939

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
1.25-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00392-1 targets PCARE in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PCARE Recombinant antibody, PBS Only (Capture) 83364-3-PBS (240352B6). 100 µg. Concentration 1 mg/mL.

Detection antibody: PCARE Recombinant antibody, PBS Only (Capture/Detector) 83364-1-PBS (240352A7). 100 µg. Concentration 1 mg/mL.

Alternative PCARE matched antibody pairs: MP00392-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00392-1, PCARE Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83364-3-PBS. Detection antibody: 83364-1-PBS. Standard: Ag34774. Range: 1.25-40 ng/mL.