

SLC6A7 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00399-2

Capture Antibody Information

Catalog Number:
83412-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240395D4
Reactivity:
human
GenBank:
BC093785
Immunogen Catalog Number:
Ag34413

Conjugate:
Unconjugated
Full name:
solute carrier family 6
(neurotransmitter transporter, L-
proline), member 7
Gene ID:
6534

Detection Antibody Information

Catalog Number:
83412-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240395E7
Reactivity:
human
GenBank:
BC093785
Immunogen Catalog Number:
Ag34413

Conjugate:
Unconjugated
Full name:
solute carrier family 6
(neurotransmitter transporter, L-
proline), member 7
Gene ID:
6534

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-80 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP00399-2 targets SLC6A7 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SLC6A7 Recombinant antibody, PBS Only (Capture) 83412-1-PBS (240395D4). 100 µg.
Concentration 1 mg/mL.

Detection antibody: SLC6A7 Recombinant antibody, PBS Only (Detector) 83412-2-PBS (240395E7). 100 µg.
Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00399-2, SLC6A7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83412-1-PBS. Detection antibody: 83412-2-PBS. Standard: Ag34413. Range: 1.25-80 ng/mL.