

Mouse IL-22 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00410-1

Capture Antibody Information

Catalog Number:

83405-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240251F9

Reactivity:

mouse

GenBank:

NM_016971

Conjugate:

Unconjugated

Full name:

interleukin 22

Gene ID:

50929

Detection Antibody Information

Catalog Number:

83405-3-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240251C8

Reactivity:

mouse

GenBank:

NM_016971

Conjugate:

Unconjugated

Full name:

interleukin 22

Gene ID:

50929

Applications

Tested Applications:

Cytometric bead array

Range:

0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00410-1 targets IL-22 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse IL22 Recombinant antibody, PBS Only (Capture/Detector) 83405-2-PBS (240251F9). 100 μ g. Concentration 1 mg/mL.Detection antibody: Mouse IL22 Recombinant antibody, PBS Only (Detector) 83405-3-PBS (240251C8). 100 μ g. Concentration 1 mg/mL.

Alternative IL-22 matched antibody pairs: MP00410-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00410-1, MOUSE IL-22 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83405-2-PBS. Detection antibody: 83405-3-PBS. Standard: Eg0405. Range: 0.313-40 ng/mL