

# Mouse IL-22 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00410-2

## Capture Antibody Information

**Catalog Number:**  
83405-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240251A12  
**Reactivity:**  
mouse  
**GenBank:**  
NM\_016971

**Conjugate:**  
Unconjugated  
**Full name:**  
interleukin 22  
**Gene ID:**  
50929

## Detection Antibody Information

**Catalog Number:**  
83405-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240251F9  
**Reactivity:**  
mouse  
**GenBank:**  
NM\_016971

**Conjugate:**  
Unconjugated  
**Full name:**  
interleukin 22  
**Gene ID:**  
50929

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00410-2 targets IL-22 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse IL22 Recombinant antibody, PBS Only (Capture) 83405-1-PBS (240251A12). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse IL22 Recombinant antibody, PBS Only (Capture/Detector) 83405-2-PBS (240251F9). 100 µg. Concentration 1 mg/mL.

Alternative IL-22 matched antibody pairs: MP00410-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00410-2, MOUSE IL-22 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83405-1-PBS. Detection antibody: 83405-2-PBS. Standard: Eg0405. Range: 0.625-40 ng/mL