

MOSPD2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00419-2

Capture Antibody Information

Catalog Number:
83398-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240282C11
Reactivity:
Human
GenBank:
BC030641
Immunogen Catalog Number:
Ag22206

Conjugate:
Unconjugated
Full name:
motile sperm domain containing 2
Gene ID:
158747

Detection Antibody Information

Catalog Number:
83398-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240282H6
Reactivity:
Human
GenBank:
BC030641
Immunogen Catalog Number:
Ag22206

Conjugate:
Unconjugated
Full name:
motile sperm domain containing 2
Gene ID:
158747

Applications

Tested Applications:
Cytometric bead array

Range:
0.78-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00419-2 targets MOSPD2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MOSPD2 Recombinant antibody, PBS Only (Capture/Detector) 83398-1-PBS (240282C11). 100 µg. Concentration 1 mg/mL.

Detection antibody: MOSPD2 Recombinant antibody, PBS Only (Detector) 83398-2-PBS (240282H6). 100 µg. Concentration 1 mg/mL.

Alternative MOSPD2 matched antibody pairs: MP00419-1, MP00419-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00419-2, MOSPD2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83398-1-PBS. Detection antibody: 83398-2-PBS. Standard: Ag22206. Range: 0.78-100 ng/mL.