

COL28A1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00423-2

Capture Antibody Information

Catalog Number:
83408-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240387B1

Reactivity:
Human

GenBank:
BC136892

Immunogen Catalog Number:
Ag21459

Conjugate:
Unconjugated

Full name:
collagen, type XXVIII, alpha 1

Gene ID:
340267

Detection Antibody Information

Catalog Number:
83408-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240387D3

Reactivity:
Human

GenBank:
BC136892

Immunogen Catalog Number:
Ag21459

Conjugate:
Unconjugated

Full name:
collagen, type XXVIII, alpha 1

Gene ID:
340267

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00423-2 targets COL28A1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: COL28A1 Recombinant antibody, PBS Only (Capture) 83408-1-PBS (240387B1). 100 µg. Concentration 1 mg/mL.

Detection antibody: COL28A1 Recombinant antibody, PBS Only (Detector) 83408-3-PBS (240387D3). 100 µg. Concentration 1 mg/mL.

Alternative COL28A1 matched antibody pairs: MP00423-1, MP00423-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00423-2, COL28A1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83408-1-PBS. Detection antibody: 83408-3-PBS. Standard: Ag21459. Range: 0.625-80 ng/mL.