

C9orf25 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00426-2**

Capture Antibody Information

Catalog Number:
83387-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240151G5
Reactivity:
Human
GenBank:
BC041009
Immunogen Catalog Number:
Ag23771

Conjugate:
Unconjugated
Full name:
chromosome 9 open reading frame 25
Gene ID:
203259

Detection Antibody Information

Catalog Number:
83387-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240151E11
Reactivity:
Human
GenBank:
BC041009
Immunogen Catalog Number:
Ag23771

Conjugate:
Unconjugated
Full name:
chromosome 9 open reading frame 25
Gene ID:
203259

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00426-2 targets C9orf25 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: C9orf25 Recombinant antibody, PBS Only (Capture) 83387-2-PBS (240151G5). 100 µg. Concentration 1 mg/mL.

Detection antibody: C9orf25 Recombinant antibody, PBS Only (Capture/Detector) 83387-1-PBS (240151E11). 100 µg. Concentration 1 mg/mL.

Alternative C9orf25 matched antibody pairs: MP00426-1, MP00426-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00426-2, C9orf25 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83387-2-PBS. Detection antibody: 83387-1-PBS. Standard: Ag23771. Range: 0.625-40 ng/mL.