

TSSK4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00436-3**

Capture Antibody Information

Catalog Number:
83401-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240373F5
Reactivity:
Human
GenBank:
BC111088
Immunogen Catalog Number:
Ag16205

Conjugate:
Unconjugated
Full name:
testis-specific serine kinase 4
Gene ID:
283629

Detection Antibody Information

Catalog Number:
83401-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240373F10
Reactivity:
Human
GenBank:
BC111088
Immunogen Catalog Number:
Ag16205

Conjugate:
Unconjugated
Full name:
testis-specific serine kinase 4
Gene ID:
283629

Applications

Tested Applications:
Cytometric bead array

Range:
6.25-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00436-3 targets TSSK4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TSSK4 Recombinant antibody, PBS Only (Capture) 83401-2-PBS (240373F5). 100 μ g. Concentration 1 mg/mL.

Detection antibody: TSSK4 Recombinant antibody, PBS Only (Detector) 83401-4-PBS (240373F10). 100 μ g. Concentration 1 mg/mL.

Alternative TSSK4 matched antibody pairs: MP00436-1, MP00436-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00436-3, TSSK4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83401-2-PBS. Detection antibody: 83401-4-PBS. Standard: Ag16205. Range: 6.25-100 ng/mL