

# FOXC2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00470-1

## Capture Antibody Information

**Catalog Number:**  
83476-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240403C6  
**Reactivity:**  
human  
**GenBank:**  
NM\_005251  
**Immunogen Catalog Number:**  
Ag33121

**Conjugate:**  
Unconjugated  
**Full name:**  
forkhead box C2 (MFH-1, mesenchyme forkhead 1)  
**Gene ID:**  
2303

## Detection Antibody Information

**Catalog Number:**  
83476-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240403B2  
**Reactivity:**  
human  
**GenBank:**  
NM\_005251  
**Immunogen Catalog Number:**  
Ag33121

**Conjugate:**  
Unconjugated  
**Full name:**  
forkhead box C2 (MFH-1, mesenchyme forkhead 1)  
**Gene ID:**  
2303

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.313-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00470-1 targets FOXC2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FOXC2 Recombinant antibody, PBS Only (Capture) 83476-3-PBS (240403C6). 100 µg. Concentration 1 mg/mL.

Detection antibody: FOXC2 Recombinant antibody, PBS Only (Capture/Detector) 83476-1-PBS (240403B2). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00470-1, FOXC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83476-3-PBS. Detection antibody: 83476-1-PBS. Standard: Ag33121. Range: 0.313-40 ng/mL