

IKZF4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00475-2

Capture Antibody Information

Catalog Number:
83471-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240357F10
Reactivity:
Human
GenBank:
BC160169
Immunogen Catalog Number:
Ag34767

Conjugate:
Unconjugated
Full name:
IKAROS family zinc finger 4 (Eos)
Gene ID:
64375

Detection Antibody Information

Catalog Number:
83471-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240357C8
Reactivity:
Human
GenBank:
BC160169
Immunogen Catalog Number:
Ag34767

Conjugate:
Unconjugated
Full name:
IKAROS family zinc finger 4 (Eos)
Gene ID:
64375

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00475-2 targets IKZF4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IKZF4 Recombinant antibody, PBS Only (Capture) 83471-4-PBS (240357F10). 100 µg. Concentration 1 mg/mL.

Detection antibody: IKZF4 Recombinant antibody, PBS Only (Detector) 83471-1-PBS (240357C8). 100 µg. Concentration 1 mg/mL.

Alternative IKZF4 matched antibody pairs: MP00475-1, MP00475-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00475-2, IKZF4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83471-4-PBS. Detection antibody: 83471-1-PBS. Standard: Ag34767. Range: 0.625-40 ng/mL