

FAM116A Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00477-2**

Capture Antibody Information

Catalog Number:
83473-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240486A7
Reactivity:
Human
GenBank:
BC040291
Immunogen Catalog Number:
Ag22052

Conjugate:
Unconjugated
Full name:
family with sequence similarity 116, member A
Gene ID:
201627

Detection Antibody Information

Catalog Number:
83473-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240486B10
Reactivity:
Human
GenBank:
BC040291
Immunogen Catalog Number:
Ag22052

Conjugate:
Unconjugated
Full name:
family with sequence similarity 116, member A
Gene ID:
201627

Applications

Tested Applications:
Cytometric bead array

Range:
0.78-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00477-2 targets FAM116A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FAM116A Recombinant antibody, PBS Only (Capture) 83473-1-PBS (240486A7). 100 µg. Concentration 1 mg/mL.

Detection antibody: FAM116A Recombinant antibody, PBS Only (Detector) 83473-2-PBS (240486B10). 100 µg. Concentration 1 mg/mL.

Alternative FAM116A matched antibody pairs: MP00477-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00477-2, FAM116A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83473-1-PBS. Detection antibody: 83473-2-PBS. Standard: Ag22052. Range: 0.78-100 ng/mL.