

ARID4A Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00489-2

Capture Antibody Information

Catalog Number: 83489-3-PBS	Clone ID: 240436C9	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: AT rich interactive domain 4A (RBP1-like)
Isotype: IgG	GenBank: NM_002892.3	Gene ID: 5926
Purification Method: Protein A purification	Immunogen Catalog Number: Ag33114	

Detection Antibody Information

Catalog Number: 83489-2-PBS	Clone ID: 240436D1	Conjugate: Unconjugated
Host: Rabbit	Reactivity: human	Full name: AT rich interactive domain 4A (RBP1-like)
Isotype: IgG	GenBank: NM_002892.3	Gene ID: 5926
Purification Method: Protein A purification	Immunogen Catalog Number: Ag33114	

Applications

Tested Applications: Cytometric bead array	Range: 1.25-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP00489-2 targets ARID4A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ARID4A Recombinant antibody, PBS Only 83489-3-PBS (240436C9). 100 µg. Concentration 1 mg/mL.

Detection antibody: ARID4A Recombinant antibody, PBS Only (Detector) 83489-2-PBS (240436D1). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00489-2, ARID4A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83489-3-PBS. Detection antibody: 83489-2-PBS. Standard: Ag33114. Range: 1.25-40 ng/mL.