

# CD86 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00511-2

## Capture Antibody Information

**Catalog Number:**

83523-3-PBS

**Host:**

Rabbit

**Isotype:**

IgG

**Purification Method:**

Protein A purification

**Clone ID:**

240429D3

**Reactivity:**

human

**GenBank:**

NM\_175862

**Conjugate:**

Unconjugated

**Full name:**

CD86 molecule

**Gene ID:**

942

## Detection Antibody Information

**Catalog Number:**

83523-2-PBS

**Host:**

Rabbit

**Isotype:**

IgG

**Purification Method:**

Protein A purification

**Clone ID:**

240429D8

**Reactivity:**

human

**GenBank:**

NM\_175862

**Conjugate:**

Unconjugated

**Full name:**

CD86 molecule

**Gene ID:**

942

## Applications

**Tested Applications:**

Cytometric bead array

**Range:**

0.313-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00511-2 targets CD86 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CD86 Recombinant antibody, PBS Only (Capture) 83523-3-PBS (240429D3). 100 µg. Concentration 1 mg/mL.

Detection antibody: CD86 Recombinant antibody, PBS Only (Detector) 83523-2-PBS (240429D8). 100 µg. Concentration 1 mg/mL.

Alternative CD86 matched antibody pairs: MP00195-1, MP00195-2, MP00195-3, MP00511-1, MP00511-3, MP50056-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage


**Storage:**

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C****Storage buffer:**

PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00511-2, CD86 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83523-3-PBS. Detection antibody: 83523-2-PBS. Standard: Eg0966. Range: 0.313-40 ng/mL.