

SYNPO2L Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00513-2**

Capture Antibody Information

Catalog Number: 83532-3-PBS	Clone ID: 240523F12	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: synaptopodin 2-like
Isotype: IgG	GenBank: BC142635	Gene ID: 79933
Purification Method: Protein A purification	Immunogen Catalog Number: Ag34027	

Detection Antibody Information

Catalog Number: 83532-4-PBS	Clone ID: 240523H7	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: synaptopodin 2-like
Isotype: IgG	GenBank: BC142635	Gene ID: 79933
Purification Method: Protein A purification	Immunogen Catalog Number: Ag34027	

Applications

Tested Applications: Cytometric bead array	Range: 1.56-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00513-2 targets SYNPO2L in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SYNPO2L Recombinant antibody, PBS Only (Capture) 83532-3-PBS (240523F12). 100 µg. Concentration 1 mg/mL.

Detection antibody: SYNPO2L Recombinant antibody, PBS Only (Detector) 83532-4-PBS (240523H7). 100 µg. Concentration 1 mg/mL.

Alternative SYNPO2L matched antibody pairs: MP00513-1, MP00513-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00513-2, SYNPO2L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83532-3-PBS. Detection antibody: 83532-4-PBS. Standard: Ag34027. Range: 1.56-100 ng/mL.