

LDHC Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00523-2**

Capture Antibody Information

Catalog Number:

83526-1-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240502B12

Reactivity:

Human

GenBank:

NM_002301

Conjugate:

Unconjugated

Full name:

lactate dehydrogenase C

Gene ID:

3948

Detection Antibody Information

Catalog Number:

83526-3-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240502E10

Reactivity:

Human

GenBank:

NM_002301

Conjugate:

Unconjugated

Full name:

lactate dehydrogenase C

Gene ID:

3948

Applications

Tested Applications:

Cytometric bead array

Range:

39-2500 pg/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00523-2 targets LDHC in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: LDHC Recombinant antibody, PBS Only (Capture) 83526-1-PBS (240502B12). 100 µg. Concentration 1 mg/mL.

Detection antibody: LDHC Recombinant antibody, PBS Only (Detector) 83526-3-PBS (240502E10). 100 µg. Concentration 1 mg/mL.

Alternative LDHC matched antibody pairs: MP00523-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:

Store at -80°C.

Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00523-2, LDHC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83526-1-PBS. Detection antibody: 83526-3-PBS. Standard: SY00609. Range: 39-2500 pg/mL.