

CD138/Syndecan-1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00527-1**

Capture Antibody Information

Catalog Number:

83551-2-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240422C10

Reactivity:

human

GenBank:

BC008765

Conjugate:

Unconjugated

Full name:

syndecan 1

Gene ID:

6382

Detection Antibody Information

Catalog Number:

83551-4-PBS

Host:

Rabbit

Isotype:

IgG

Purification Method:

Protein A purification

Clone ID:

240422C12

Reactivity:

human

GenBank:

BC008765

Conjugate:

Unconjugated

Full name:

syndecan 1

Gene ID:

6382

Applications

Tested Applications:

Cytometric bead array

Range:

0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00527-1 targets CD138/Syndecan-1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CD138/Syndecan-1 Recombinant antibody, PBS Only (Capture/Detector) 83551-2-PBS (240422C10). 100 µg. Concentration 1 mg/mL.

Detection antibody: CD138/Syndecan-1 Recombinant antibody, PBS Only (Detector) 83551-4-PBS (240422C12). 100 µg. Concentration 1 mg/mL.

Alternative CD138/Syndecan-1 matched antibody pairs: MP00527-2, MP00527-3, MP00990-1, MP00990-2, MP00990-3, MP50044-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00527-1, CD138/Syndecan-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83551-2-PBS. Detection antibody: 83551-4-PBS. Standard: Eg0132. Range: 0.313-40 ng/mL