

Mouse Chemerin Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00537-1**

Capture Antibody Information

Catalog Number: 83550-2-PBS	Clone ID: 240355D4	Conjugate: Unconjugated
Host: Rabbit	Reactivity: mouse	Full name: retinoic acid receptor responder (tazarotene induced) 2
Isotype: IgG	GenBank: NM_001347168.1	Gene ID: 71660
Purification Method: Protein A purification	Immunogen Catalog Number: Eg0875	

Detection Antibody Information

Catalog Number: 83550-1-PBS	Clone ID: 240355B3	Conjugate: Unconjugated
Host: Rabbit	Reactivity: mouse	Full name: retinoic acid receptor responder (tazarotene induced) 2
Isotype: IgG	GenBank: NM_001347168.1	Gene ID: 71660
Purification Method: Protein A purification	Immunogen Catalog Number: Eg0875	

Applications

Tested Applications: Cytometric bead array	Range: 0.313-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00537-1 targets Chemerin in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse Chemerin Recombinant antibody, PBS Only (Detector) 83550-2-PBS (240355D4). 100 µg. Concentration 1 mg/mL.

Detection antibody: Mouse Chemerin Recombinant antibody, PBS Only (Capture) 83550-1-PBS (240355B3). 100 µg. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00537-1, MOUSE Chemerin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83550-2-PBS. Detection antibody: 83550-1-PBS. Standard: Eg0875. Range: 0.313-40 ng/mL