

ARID3B Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00548-2**

Capture Antibody Information

Catalog Number:
83560-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240466F5

Reactivity:
Human

GenBank:
BC041792

Immunogen Catalog Number:
Ag10545

Conjugate:
Unconjugated

Full name:
AT rich interactive domain 3B (BRIGHT-like)

Gene ID:
10620

Detection Antibody Information

Catalog Number:
83560-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240466F1

Reactivity:
Human

GenBank:
BC041792

Immunogen Catalog Number:
Ag10545

Conjugate:
Unconjugated

Full name:
AT rich interactive domain 3B (BRIGHT-like)

Gene ID:
10620

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00548-2 targets ARID3B in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ARID3B Recombinant antibody, PBS Only (Capture/Detector) 83560-2-PBS (240466F5). 100 μ g. Concentration 1 mg/mL.

Detection antibody: ARID3B Recombinant antibody, PBS Only (Detector) 83560-3-PBS (240466F1). 100 μ g. Concentration 1 mg/mL.

Alternative ARID3B matched antibody pairs: MP00548-1, MP00548-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00548-2, ARID3B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83560-2-PBS. Detection antibody: 83560-3-PBS. Standard: Ag10545. Range: 0.781-100 ng/mL.