

DHRS3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00554-2**

Capture Antibody Information

Catalog Number:
83581-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240545F6
Reactivity:
Human
GenBank:
BC002730
Immunogen Catalog Number:
Ag7589

Conjugate:
Unconjugated
Full name:
dehydrogenase/reductase (SDR family) member 3
Gene ID:
9249

Detection Antibody Information

Catalog Number:
83581-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240545C4
Reactivity:
Human
GenBank:
BC002730
Immunogen Catalog Number:
Ag7589

Conjugate:
Unconjugated
Full name:
dehydrogenase/reductase (SDR family) member 3
Gene ID:
9249

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00554-2 targets DHRS3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DHRS3 Recombinant antibody, PBS Only (Capture) 83581-2-PBS (240545F6). 100 µg. Concentration 1 mg/mL.

Detection antibody: DHRS3 Recombinant antibody, PBS Only (Capture/Detector) 83581-1-PBS (240545C4). 100 µg. Concentration 1 mg/mL.

Alternative DHRS3 matched antibody pairs: MP00554-1, MP00554-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00554-2, DHR53 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83581-2-PBS. Detection antibody: 83581-1-PBS. Standard: Ag7589. Range: 0.313-40 ng/mL