

# MLL5 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00627-2**

## Capture Antibody Information

**Catalog Number:**  
83661-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240664A9  
**Reactivity:**  
Human  
**GenBank:**  
BC062583  
**Immunogen Catalog Number:**  
Ag6044

**Conjugate:**  
Unconjugated  
**Full name:**  
myeloid/lymphoid or mixed-lineage leukemia 5 (trithorax homolog, Drosophila)  
**Gene ID:**  
55904

## Detection Antibody Information

**Catalog Number:**  
83661-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240664A6  
**Reactivity:**  
Human  
**GenBank:**  
BC062583  
**Immunogen Catalog Number:**  
Ag6044

**Conjugate:**  
Unconjugated  
**Full name:**  
myeloid/lymphoid or mixed-lineage leukemia 5 (trithorax homolog, Drosophila)  
**Gene ID:**  
55904

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-80 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00627-2 targets MLL5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MLL5 Recombinant antibody, PBS Only (Capture) 83661-3-PBS (240664A9). 100 µg. Concentration 1 mg/mL.

Detection antibody: MLL5 Recombinant antibody, PBS Only (Capture/Detector) 83661-1-PBS (240664A6). 100 µg. Concentration 1 mg/mL.

Alternative MLL5 matched antibody pairs: MP00627-1, MP00627-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

# Selected Validation Data



Cytometric bead array standard curve of MP00627-2, MLL5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83661-3-PBS. Detection antibody: 83661-1-PBS. Standard: Ag6044. Range: 0.625-80 ng/mL