

GBP4 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00655-2

Capture Antibody Information

Catalog Number:
83699-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240730F11
Reactivity:
Human
GenBank:
BC070055
Immunogen Catalog Number:
Ag12006

Conjugate:
Unconjugated
Full name:
guanylate binding protein 4
Gene ID:
115361

Detection Antibody Information

Catalog Number:
83699-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240730A12
Reactivity:
Human
GenBank:
BC070055
Immunogen Catalog Number:
Ag12006

Conjugate:
Unconjugated
Full name:
guanylate binding protein 4
Gene ID:
115361

Applications

Tested Applications:
Cytometric bead array

Range:
3.13-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00655-2 targets GBP4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GBP4 Recombinant antibody, PBS Only (Capture) 83699-2-PBS (240730F11). 100 µg. Concentration 1 mg/mL.

Detection antibody: GBP4 Recombinant antibody, PBS Only (Detector) 83699-3-PBS (240730A12). 100 µg. Concentration 1 mg/mL.

Alternative GBP4 matched antibody pairs: MP00655-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00655-2, GBP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83699-2-PBS. Detection antibody: 83699-3-PBS. Standard: Ag12006. Range: 3.13-100 ng/mL