

MTMR9 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: **MP00672-3**

Capture Antibody Information

Catalog Number:
83694-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240167E6
Reactivity:
Human
GenBank:
BC022003
Immunogen Catalog Number:
Ag34134

Conjugate:
Unconjugated
Full name:
myotubularin related protein 9
Gene ID:
66036

Detection Antibody Information

Catalog Number:
83694-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240167E10
Reactivity:
Human
GenBank:
BC022003
Immunogen Catalog Number:
Ag34134

Conjugate:
Unconjugated
Full name:
myotubularin related protein 9
Gene ID:
66036

Applications

Tested Applications:
Cytometric bead array

Range:
0.78-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00672-3 targets MTMR9 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MTMR9 Recombinant antibody, PBS Only (Capture/Detector) 83694-1-PBS (240167E6). 100 µg. Concentration 1 mg/mL.

Detection antibody: MTMR9 Recombinant antibody, PBS Only (Detector) 83694-3-PBS (240167E10). 100 µg. Concentration 1 mg/mL.

Alternative MTMR9 matched antibody pairs: MP00672-1, MP00672-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00672-3, MTMR9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83694-1-PBS. Detection antibody: 83694-3-PBS. Standard: Ag34134 Range: 0.78-100 ng/mL